

IOWA STATE UNIVERSITY

AIRFOIL WAKE MEASUREMENTS AND
CALIBRATION OF A HOT WIRE ANEMOMETER
PRE-LABORATORY

AER E 344 - PRE-LAB 06 - AIRFOIL WAKE MEASUREMENTS AND
CALIBRATION OF A HOT WIRE ANEMOMETER

SECTION 3 GROUP 3

MATTHEW MEHRTENS
JACK MENDOZA
KYLE OSTENDORF
GABRIEL PEDERSON
LUCAS TAVARES VASCONCELLOS
DREW TAYLOR

PROFESSOR

HUI HU, PhD

College of Engineering
Aerospace Engineering
Aerodynamics and Propulsion Laboratory

AMES, FEBRUARY 2024

ANSWERS

1.1 Question 1

You should review and understand the fundamental technical basis of thermal-based anemometry techniques.

Thermal anemometers use two thin wire thermometers to measure the flow velocity of a gas. One wire measures the temperature of the fluid as a reference while the other wire is heated using an electrical current. The flow removes heat from the heated wire and causes a cooling effect which increases as the flow increases. The temperatures of the wires are not measured directly, but are calculated using the relationship with resistance.

1.2 Question 2

You should understand the differences between CCA and CTA approaches in using hotwire anemometry system for flow velocity measurements.

CCA - In constant-current anemometry, the current passing through the heated wire remains constant. The fluid flow fluctuations cause change in the wire temperature which can be used to calculate the flow velocity in a first order differential equation.

CTA – Constant-temperature anemometry keeps the hot wire the same temperature using an electric feedback system to adjust the voltage going to the hot wire in response to the cooling effect of fluid flow.

1.3 Procedures

You should review the recorded video of the AER E 344 Pre-lab 06.

1.3.1 Airfoil Wake Measurements

1. Calibrate the instruments when the wind tunnel is at 0 Hz
2. Set the velocity at 10 Hz to 15 Hz. Wait for the velocity to settle.
3. Set the angle of attack (AoA) to -4° .
4. Move the rake to cover the entire wake of the airfoil as necessary
5. Acquire and save the data (it takes approximately 5 s)
6. Repeat **Step 3** through **Step 5** using the following AoAs: 0° , 4° , 6° , 8° , 10° , 12° and 16°

1.3.2 Hot-Wire Anemometer Calibration

1. Get the voltage data when the velocity is 0 m/s
2. Set the wind tunnel to 5 Hz and record the data. This includes the voltage data given by the computer and the pressure data from the Mensor manometer.
3. Repeat **Step 2** over a range of frequencies from 5 Hz to 35 Hz, incrementing by 5 Hz each time
4. Approximate the data to a 4th degree function